

Shimadzu EcoCells

EcoCells from Shimadzu is a economical range of Cuvette cells that offers different light path ranging from 10mm to 100mm. These cuvette cells are applicable on every spectrophotometer and photometer instrument in the market. These quartz cuvettes includes PTFE lid.

Approximately 80% transmission on empty cell for EcoCells at 200nm.



Light Path/mm	External Dimension	Single cell	Matched pair
10	45 x 12.5 x 12.5	226-87010-91	226-87010-92
40	45 x 12.5 x 42.5	226-87040-91	-
50	45 x 12.5 x 52.5	226-87050-91	226-87050-92
100	45 x 12.5 x 102.5	-	226-87100-92